## **Electro-Optic Coefficients Calculation** (PyEOC)

Implementation of the method published by Cuniot-Ponsard et al. in JAP 109, 014107 (2011) <a href="https://doi.org/10.1063/1.3514083">https://doi.org/10.1063/1.3514083</a>

Pr. Sidi Hamady Université de Lorraine, France

sidi.hamady@univ-lorraine.fr

Released under the MIT license ( $\underline{\text{https://opensource.org/licenses/MIT}}$ ) See Copyright Notice in COPYRIGHT

Reference: Sidi Ould Saad Hamady, "PyEOC: a Python Code for Determining Electro-Optic Coefficients of Thin-Film Materials", 2022.

## **Presentation and Requirements**

PyEOC calculates the electro-optic coefficients of an optical material by implementing the method published by Cuniot-Ponsard et al. in JAP 109, 014107 (2011) <a href="https://doi.org/10.1063/1.3514083">https://doi.org/10.1063/1.3514083</a>.





Figure 1. The standard multilayer structure used to measure the electro-optic (and converse piezoelectric) coefficients of a thin-film material. The film to characterize is sandwiched between two thin metallic contacts (e.g. gold, platinum): a semitransparent top contact and a thick bottom contact. The voltage is applied between these two contacts and the reflectance measured with respect to the incident angle 0.



**Figure 2.** Intensity variation with the position in the Pt/SBN/Pt structure in both TE and TM polarization. The figures on the right represent a zoom in the active layer.




Figure 3. Experimental and theoretical calculated static reflectivity R (top-left), dynamic reflectivity  $\Delta R$  (bottom-left) and the calculated TE derivatives (right). The symbols indicate the experimental data curves and the lines indicate the fitted curves.

## To install PyEOC:

just download it from github: <a href="https://github.com/sidihamady/PyEOC">https://github.com/sidihamady/PyEOC</a>

unzip and use.

You can also install the PyEOC module to use as usually done in Python:

```
python3 setup.py install
```

The distribution mainly includes:

Three main Python files:

- PyEOC.py implementing the program core functionality in the module classe.
- <a href="mailto:tmmCore.py">tmmCore.py</a> S J Byrnes' transfer-matrix-method code

  (<a href="https://arxiv.org/abs/1603.02720">https://arxiv.org/abs/1603.02720</a>)
- <u>Test.py</u> an example with a SBN structure by Cuniot-Ponsard et al. in JAP 109, 014107 (2011)

The basic requirements are found in any Linux distribution (and easily installed for Windows):

- Python version 2.7.x or later
- numpy version 1.5 or later
- scipy version 0.13.1 or later
- matplotlib version 1.3.x or later

PS: for Windows, you can download a complete Python distribution from <a href="https://www.anaconda.com/distribution/">https://www.anaconda.com/distribution/</a>

## **HowTo**

You can use the included Test.py and adapt it to your needs:

```
# -*- coding: utf-8 -*-
import sys, os
#sys.path.insert(0, "/path/to/PyEOC")
from PyEOC import * # import PyEOC core class
Structure = PyEOC(
    'SBN', # structure name included in the PyEOC class
   # measurement data: static reflectivity vs angle (TE and TM)
                       dynamic reflectivity vs angle (TE and TM)
                       four files with tab-separated columns
   # 'SBN' data extracted from Cuniot-Ponsard et al. in JAP 109, 014107 (2011)
   'SBN_Reflectivity_TE.txt',
                                 'SBN_Reflectivity_TM.txt',
    'SBN_Reflectivity_Dyn_TE.txt', 'SBN_Reflectivity_Dyn_TM.txt'
Structure.wavelength = 633 # laser wavelength in nm
Structure.voltage = 1.0
                          # applied voltage amplitude in volts
# the incident angle theta starting three values and range
# the choosen theta values should correspond to a "smooth" and "different" part of the
DR and delta_R / delta_? data
Structure.theta_manual = [35.0 * Structure.toradian, 40.0 * Structure.toradian, 45.0
* Structure.toradian]
Structure.thetaDelta = 2.0 * Structure.toradian
```

```
Structure.thetaStart = 30.0 * Structure.toradian

Structure.thetaEnd = 70.0 * Structure.toradian

Structure.thickness[1] = 22.6 # Pt thickness (nm)

Structure.thickness[2] = 758 # SBN thickness (nm)

Structure.refractiveindexo[2] = 2.30 + 0.0515624j

Structure.refractiveindexe[2] = 2.26 + 0.0515624j

Structure.fit_dynamic = True # fit the dynamic reflectivity?

Structure.fit(report = True) # start fitting and report

Structure.plot() # plot the fitted curves

#Structure.plotpoynting() # plot the intensity profile
```